



DISPOSABLE  
LAPAROSCOPY



Founded over 20 years ago, **Neomedic International** is a company entirely focused on providing the medical field with added value surgical methods, solutions and innovative products, which improve physician quality of work as well as patient quality of life.

The greatest value for **Neomedic International** are the people. Customers and patients are at the forefront of every decision and action we take; from product design to long-term product performance, clinicians and patients always come first.

**Neomedic International** works alongside medical professionals who are in the OR on a daily basis, to further comprehend of the needs of surgeons but most importantly, to better treat patient conditions. Our passion is to listen to our customers and help transform their needs into useful products for more successful patient care.








**Neomedic International** has a fully integrated manufacturing facility with technologies that also act as a flexible R&D resource.

Words like responsibility, efficiency, and innovation are all part of **Neomedic International's** philosophy. We are proud to be your partner.





## NEOLAP LAPAROSCOPY

### DISPOSABLE LAPAROSCOPY

	<p><b>TROCARS</b></p> <ul style="list-style-type: none"> <li>• Optical and bladed trocars</li> <li>• Balloon trocars</li> </ul>
	<p><b>MONOPOLAR INSTRUMENTS</b></p> <ul style="list-style-type: none"> <li>• Monopolar instruments 3mm Ø, standard and pediatric length</li> <li>• Monopolar instruments 5mm Ø, standard and bariatric length</li> </ul>
	<p><b>BIPOLAR INSTRUMENTS</b></p> <ul style="list-style-type: none"> <li>• Bipolar instruments 3mm Ø and 5mm Ø</li> </ul>
	<p><b>CLIP APPLIERS AND POLYMER LIGATION CLIPS</b></p> <ul style="list-style-type: none"> <li>• Clip applicators</li> <li>• Polymer clips</li> </ul>
	<p><b>SUCTION/ IRRIGATION SYSTEMS</b></p> <ul style="list-style-type: none"> <li>• Disposable systems, ergonomic and trumpet</li> <li>• Monopolar electrodes and cannulas</li> </ul>
	<p><b>RETRIEVAL BAGS</b></p>
	<p><b>VERESS NEEDLES</b></p>

### OTHER NEOLAP RANGE PRODUCTS

	<p><b>UTERUS MANIPULATOR</b></p>
	<p><b>CLIP APPLIERS FOR OPEN SURGERY</b> <b>ELICK EVACUATORS</b></p>

# OPTICAL TROCARS

3MM, 5MM 10MM AND 12MM DIAMETER TROCARS

## OPTICAL TROCARS



### 3mm Ø

101Y.510– 3mm. Ø trocar with dilating tip. Purple cannula. 15u/box.

101Y.510B– 3mm. Ø purple cannula. 15u/box

### 5mm Ø

101Y.511– 5mm Ø trocar with dilating tip. Purple cannula. 15u/box

101Y.511B– 5mm purple cannula. 15u/box

101Y.517– KIT: 1unit 5mm Ø trocar + 2 units purple cannulas. 10u/box

### 10mm Ø

101Y.512– 10mm Ø trocar with dilating tip. Purple cannula. 15u/box

101Y.512B– 10mm Ø purple cannula. 15u/box

101Y.518– KIT: 1unit 10mm Ø trocar + 2 units purple cannulas. 10u/box

### 12mm Ø

101Y.513– 12mm Ø trocar with dilating tip. Purple cannula. 15u/box

101Y.513B– 10mm Ø purple cannula. 15u/box

101Y.519– KIT: 1unit 10mm Ø trocar + 2 units purple cannulas. 10u/box

# BLADED TROCARS

5MM, 10MM AND 12MM DIAMETER TROCARS

## BLADED TROCARS

DETACHABLE

TRANSPARENT CANNULA

UNIVERSAL SEALING

No reducer needed

THREADED DESIGN

Maximizes abdominal  
wall retention

AUTORSHIELD BLADE

Safe insertion

### Ø 5mm

101Y.524— 5mm Ø trocar with autoshield tip. Purple cannula. 15u/box

101Y.511B— 5mm. Ø purple cannula. 15u/box

101Y.421 KIT: 1unit 5mm Ø trocar + 2 units purple cannulas. 10u/box

### Ø 10mm

101Y.522— 10mm Ø trocar with autoshield tip. Purple cannula. 15u/box

101Y.512B— 10mm. Ø purple cannula. 15u/box

101Y.422— KIT: 1unit 10mm Ø trocar + 2 units purple cannulas. 10u/box

### Ø 12mm

101Y.523— 12mm Ø trocar with autoshield tip.

Purple cannula. 15u/box

101Y.513B— 12mm. Ø purple cannula. 15u/box

101Y.423— KIT: 1unit 12mm Ø trocar + 2 units purple  
cannulas. 10u/box

# BALLOON OPTICAL TROCARS

5MM, 10MM AND 12MM DIAMETER TROCARS

DETACHABLE

TRANSPARENT CANNULA

SMOOTH SURFACE



UNIVERSAL SEALING  
No reducer needed




FIXATION BALLON

Transparent to  
maximize visibility



RETENTION DISK

Sliding piece, fixes  
perfectly providing pressure  
on the incision site



## Ø 5mm

BOTY.005— 5mm Ø trocar with dilating tip. Balloon volume: 1,5ml. 8 u/box

## Ø 10mm

BOTY.010— 10mm Ø trocar with dilating tip. Balloon volume: 4,0ml. 6 u/box

## Ø 12mm


BT.012— 12mm Ø trocar with dilating tip. Balloon volume: 5,0ml. 6 u/box



# MONOPOLAR AND BIPOLAR INSTRUMENTS

## STANDARD AND PEDIATRIC LENGTH LAPAROSCOPIC INSTRUMENTS

### 3MM Ø MONOPOLAR INSTRUMENTS



#### 3mm Ø– 330mm. length

- 3N.701— Curved dissector Maryland, grey color code. 14mm stainless steel tip. 12u/box
- 3N.802— Metzenbaum curved scissors, orange color code, 12mm stainless steel tip. 12u/box
- 3N.703\*— Atraumatic fenestrated grasper, blue color code. 13mm stainless steel tip. **With ratchet.** 12u/box
- 3N.704\*— Clinch grasper, yellow color code. 16mm stainless steel tip. **With ratchet.** 12u/box

#### 3mm Ø – 240mm. length (pediatric)


- 3N.701-24— Curved dissector Maryland, grey color code. 14mm stainless steel tip. 12u/box
- 3N.802-24— Metzenbaum curved scissors, orange color code, 12mm stainless steel tip. 12u/box
- 3N.703-24\*— Atraumatic fenestrated grasper, blue color code. 13mm stainless steel tip. **With ratchet.** 12u/box
- 3N.704-24\*— Clinch grasper, yellow color code. 16mm stainless steel tip. **With ratchet.** 12u/box



# MONOPOLAR AND BIPOLAR INSTRUMENTS

## STANDARD LENGTH LAPAROSCOPIC INSTRUMENTS

### 5MM Ø MONOPOLAR INSTRUMENTS



### Ø 5mm– 330mm. Length, 60° max. opening

- 101Y.701– Curved dissector Maryland, grey color code, 19mm stainless steel tip. 12u/box.
- 101Y.702– Curved Metzenbaum scissors. Black color code. 16,5mm stainless steel tip. 12u/box
- 106Y.103– Isolated Curved Metzenbaum Neolap mini-scissors. White color code. 12mm stainless steel tip. 12u/box.
- 101Y.802– Isolated Curved Metzenbaum scissors. Orange color code. 16,5mm stainless steel tip. 12u/box.
- 101Y.703\*– Atraumatic fenestrated grasper. Pale blue color code. 19mm stainless steel tip. **With ratchet.** 12u/box.
- 101Y.704\*– Clinch grasper, yellow color code. 19mm stainless steel tip. **With ratchet.** 12u/box
- 101Y.705\*– Allis grasper. 19mm stainless steel tip. **With ratchet.** 12u/box
- 101Y.706\*– Babcock grasper. 19mm stainless steel tip. **With ratchet.** 12u/box
- 101Y.707\*– Blunt grasper. 19mm stainless steel tip. **With ratchet.** 12u/box






# MONOPOLAR AND BIPOLAR INSTRUMENTS

## STANDARD LENGTH LAPAROSCOPIC INSTRUMENTS

### 5MM Ø MONOPOLAR INSTRUMENTS



### 5mm.Ø – 330mm. length—60° max. opening

**5N.701**—Curved Maryland dissector. Orange color code. 19mm stainless steel tip.12u/box

**5N.802**— Curved Metzenbaum scissors. Red color code. 12u/box

**5N.703\***— Atraumatic fenestrated double action grasper. Navy blue color code. **With ratchet.** 12u/box

**5N.704\***— Clinch grasper. Green color code. **With ratchet.** 12u/box



5N.802



5N.701



5N.703




5N.704

# MONOPOLAR AND BIPOLAR INSTRUMENTS

## BARIATRIC LENGTH LAPAROSCOPIC INSTRUMENTS

### 5MM Ø MONOPOLAR INSTRUMENTS



### 5mm Ø– 450mm. length (bariatric) 60° max. opening

101Y.701L—Curved dissector Maryland, grey color code, 19mm stainless steel tip. 12u/box.

101Y.702L— Curved Metzenbaum scissors. Black color code. 16,5mm stainless steel tip. 12u/box.

101Y.802L— Isolated Curved Metzenbaum scissors. Orange color code. 16,5mm stainless steel tip. 12u/box.

106Y.103L— Isolated Curved Metzenbaum mini-scissors. White color code. 12mm stainless steel tip. 12u/box.

101Y.703L\*— Atraumatic fenestrated grasper. Light blue color code. 19mm stainless steel tip. **With ratchet.** 12u/box.

101Y.704L\*— Clinch grasper, yellow color code. 19mm stainless steel tip. **With ratchet.** 12u/box.

101Y.705L\*— Allis grasper. 19mm stainless steel tip. **With ratchet.** 12u/box.

101Y.706L\*— Babcock grasper. 19mm stainless steel tip. **With ratchet.** 12u/box.


101Y.707L\*— Blunt grasper. 19mm stainless steel tip. **With ratchet.** 12u/box.



# MONOPOLAR AND BIPOLAR INSTRUMENTS

## STANDARD LENGTH DISPOSABLE LAPAROSCOPIC INSTRUMENTS


### 3MM AND 5MM Ø BIPOLAR INSTRUMENTS



#### 3mm Ø – 330mm length- 60° max. opening

3N.701BIPO— Maryland curved dissector, 14mm stainless steel tip. 12u/box

3N.703BIPO— Atraumatic, fenestrated grasper. 13mm stainless steel tip. 12u/box



3N.703BIPO



3N.701BIPO

#### 5mm Ø – 330mm length

109Y.401\*— Cutting and coagulating bipolar instrument. Ergonomic handle **with ratchet**.. 16,5mm stainless steel tip.  
Cable included. 1u

109Y.202\*—Maryland curved dissector. Ergonomic handle **with ratchet**. 19mm stainless steel tip. Cable included. 5u/box

109Y.201\*— Atraumatic, fenestrated grasper. Ergonomic handle **with ratchet**. 19mm stainless steel tip. Cable included. 5u/box



109Y.201



109Y.202

# REUSABLE CLIP APPLIERS


## REUSABLE CLIP APPLIERS AND DISPOSABLE POLYMER

### 5MM AND 10MM POLYMER CLIP APPLIERS



#### NON ABSORBABLE CLIPS

To ligate from 2mm. to 16mm.  
of vessel and/or tissue



#### CODE COLOR

Reusable appliers.  
Each color corresponds to  
the color of the cartridges



#### NON ABSORBABLE CLIPS

Non conductive and CT/MRI compatibles.  
For both laparoscopic and open surgery




### 5mmØ reusable appliers– 330mm length

101.113B– Laparoscopic applicator. Reusable. Size M (Green). 1u

### 10mmØ reusable appliers– 330mm length

101.113A–Laparoscopic applicator. Reusable. Size L (Purple). 1u

101.113– Laparoscopic applicator. Reusable. Size XL (Gold). 1u



### Disposable Polymer Clips

101Y.623– Polymer Clip. Size M (Green). Ligation range 3mm-10mm. 6 clips x cartridge, 14 cartridges x box.


101Y.622–Polymer Clip. Size L (Purple). Ligation range 5mm-13mm. 6 clips x cartridge, 14 cartridges x box.

101Y.621–Polymer Clip. Size XL (Gold). Ligation range 7mm-16mm. 6 clips x cartridge, 14 cartridges x box.

# SUCTION-IRRIGATION SYSTEMS

## ERGONOMIC AND TRUMPET SYSTEMS

### SUCTION-IRRIGATION SYSTEMS



**ERGONOMIC HANDPIECE**

Easy to use with both hands




**TRUMPET HANDPIECE**

Transparent handpiece for better obstructions visualization



**REUSABLE CANNULAS**

3mm, 5mm and 10mm. Ø cannulas



### Suction / Irrigation

SI07-360E— ERGONOMIC system for laparoscopic surgery. 5mm Ø and 360mm length. 12u/box

SI06-360T— TRUMPET system for laparoscopic surgery. 5mm Ø and 360mm length. 12u/box

### 3mmØ, 5mmØ and 10mmØ reusable compatible cannulas

CA03-360—3mm Ø and 330mm length stainless steel reusable cannula for minilaparoscopy surgery. 1u

CA10-360—10mmØ and 360mm length stainless steel reusable cannula for laparoscopy surgery.1u

CA05-450—5mmØ and 450mm length stainless steel reusable cannula for bariatric surgery.1u



### 3mmØ and 5mmØ disposable monopolar electrodes

LH07-330—5mm Ø suction-irrigation cannula with retractable L-Hook tip electrode for cauterization. 1u

SP07-330- 5mm Ø suction-irrigation cannula with retractable Spatula tip electrode for cauterization. 1u

3N.LHOOK- 3mm Ø and 330mm length suction-irrigation cannula with retractable L-Hook tip electrode for cauterization. 12u/box

3N.LHOOK-24- 3mm Ø and 240mm length suction-irrigation cannula with retractable L-Hook tip electrode for cauterization.12u/box

LHOOK- 5mm Ø suction-irrigation L-Hook tip electrode for cauterization. 50u /box

SH00K- 5mm Ø suction-irrigation cannula with retractable Spatula tip electrode for cauterization. 50u/ box.



# DISPOSABLE RETRIEVAL BAGS

10MM Ø RETRIEVAL BAGS AND 330MM ROD LENGTH

## AUTOMATIC RETRIEVAL BAGS



### 10 mm Ø automatic single use bags and 330mm. shaft.

**RB10.250**—Automatic retrieval bag for organ and tissue recollection. **250ml.** capacity.

Bag dimensions: 119mm x 244mm. 20u/box.

**RB10.350**—Automatic retrieval bag for organ and tissue recollection. **350ml.** capacity.

Bag dimensions: 163 x 260mm. 20u/box.

**104Y.304**—Automatic retrieval bag for organ and tissue recollection. **500ml.** capacity.

Bag dimensions: 80 x 200mm. 12u/box

**RB10.1000**—Automatic retrieval bag for organ and tissue recollection. **1000ml.** capacity.

Bag dimensions: 224 x 296mm. 20u/box.

**RB10.1200**—Automatic retrieval bag for organ and tissue recollection. **2100ml.** capacity.

Bag dimensions: 246 x 166mm. 20u/box.

# VERESS NEEDLE


VERESS NEEDLES 120MM AND 150MM LENGTH

## DISPOSABLE VERESS NEEDLES




**ERGONOMIC PC HANDLE AND VALVE**

Made of transparent TPU




**STAINLESS STEEL TROCAR**

2.1mm diameter Trocar



**RETRACTIL BLADE TIP**

Prevents damaging intraabdominal structures



## Veress insufflation needle


**101Y.550-** Veress insufflation needle, 120mm length. 10u/box

**101Y.551-** Veress insufflation needle, 150mm length. 10u/box

# UTERUS MANIPULATORS

## DOUBLE CUP DISPOSABLE UTERUS MANIPULATORS

### DISPOSABLE UTERUS MANIPULATORS



#### DOUBLE CUP DESIGN

The forward cup fits over the cervix, the back cup prevents loss of pneumo.




#### ERGONOMIC HANDLE

13cm. length handle



#### DYE TEST PORT

Dye test port for contrast injection



### Disposable 13 cm. handle uterus manipularors

**L-UM.5S**— Uterus manipulator size S. Anterior cup 29mm internal diameter and 32 mm external diameter. 15mm tall. 8u/box.

**L-UM.5M**—Uterus manipulator size M. Anterior cup 30mm internal diameter and 35 mm external diameter. 20mm tall. 8u/box.

**L-UM.5L**—Uterus manipulator size L. Anterior cup 34mm internal diameter and 37 mm external diameter. 24mm tall. 8u/box.

**L-UM.5XL**—Uterus manipulator size XL. Anterior cup 36mm internal diameter and 40 mm external diameter. 30mm tall. 8u/box.



# INSTRUMENTS

## OPEN SURGERY


### REUSABLE 5MM AND 10MM POLYMER CLIP APPLIERS

REUSABLE CLIP APPLIERS FOR OPEN SURGERY AND POLYMER CLIPS



#### NON ABSORBABLE CLIPS

To ligate from 2mm. to 16mm.  
of vessel and/or tissue



#### CODE COLOR

Reusable applicers.  
Each color corresponds to  
the color of the cartridges



#### NON ABSORBABLE CLIPS

Non conductive and CT/MRI compatibles.  
For both laparoscopic and open surgery



### 5mm. Ø Appliers

105.11B— Clip applicator for **open surgery**, 5mm Ø and **250mm length**. Reusable. size M (green). 1 u/box.

### 10mm. Ø Appliers

105.113A— Clip applicator for **open surgery**, 10mm Ø and **250mm length**. Reusable. size L (purple). 1u.

105.113— Clip applicator for **open surgery**, 10mm Ø and **250mm length**. Reusable. size XL (gold). 1u.

## ELLICK EVACUATOR

ELLICK EVACUATOR FOR PROSTATIC SURGERY



To use during prostatic surgery. For suction and irrigation of tissue, blood or stones from the urinary bladder.

280ml of capacity, sizes: 62,7mm Ø and 216,74mm.

EE.STD— Ellick Evacuator for prostatic surgery. 6u/box



Contact: [Info@neomedic.com](mailto:Info@neomedic.com)





**surgical innovation**

© 2019 Neomedic International. All rights reserved.



[www.neomedic.com](http://www.neomedic.com)