

# DISPOSABLE COVER FOR HYSTEROSCOPY

Gynko<sup>®</sup>

3MM

with 5 FR CHANNEL



## ADVANTAGES



*It guarantees the complete sterility of camera and optics*



*It allows to avoid any type of wear and tear resulting from the sterilization cycles*



*It reduces the duration of the examination procedure*



### GFRGYN300

GYNKO® for rigid optics, 2,9mm, 300mm, 0° rigid lens with Storz-like connection and camera cover.

**Box: 10 pieces**

### GFRGYN300W

Without camera cover

**Box: 10 pieces**

### GFRGYN330

GYNKO® for rigid optics, 2,9mm, 300mm, 30° rigid lens with Storz-like connection and camera cover.

**Box: 10 pieces**

### GFRGYN330W

Without camera cover

**Box: 10 pieces**

### GFRGYN430

GYNKO® for rigid optics, 4mm, 300mm, 0°- 30° smooth lens with Storz-like connection and camera cover.

**Box: 10 pieces**

### GFRGYN430W

Without camera cover

**Box: 10 pieces**

CND: T030102

RDM: 641108

Distributed by:



# Gynko®